

Product Datasheet

IgG antibody (HRP) (orb436473)

Description	Donkey polyclonal antibody against IgG conjugated to HRP
Species/Host	Donkey
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Target	IgG
Preservatives	1.5% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Concentration	Current, batch-specific concentration 0.8 mg/ml after reconstitution
Storage	-20°C after reconstitution
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.