

Product Datasheet

IgG1 antibody (HRP) (orb436475)

Description	Rabbit polyclonal antibody against IgG1 conjugated to HRP
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	DOT, ELISA, IHC-Fr, WB
Immunogen	Mouse IgG1 gamma 1 heavy chain
Target	IgG1
Preservatives	0.01% gentamicin sulphate1% Bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.