

Product Datasheet

IgM antibody (FITC) (orb436479)

Description	Goat polyclonal antibody against IgM conjugated to FITC
Species/Host	Goat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Target	IgM
Preservatives	0.05% Sodium Azide (NaN ₃)1.5% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - lyophilized
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.