

Product Datasheet

IgG antibody (orb436482)

Description	Goat polyclonal antibody to IgG
Species/Host	Goat
Reactivity	Hamster
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Hamster IgG.
Target	IgG
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	F(ab') ₂ fragment of purified IgG - liquid
Concentration	0.8 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.