

Product Datasheet

IgG antibody (AP) (orb436528)

Description	Goat polyclonal antibody against IgG conjugated to AP
Species/Host	Goat
Reactivity	Human
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	Human IgG F(ab') ₂
Target	IgG
Preservatives	0.1% Sodium Azide (NaN ₃)1.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified Ig conjugated to Alkaline Phosphatase - liquid
Concentration	0.6mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	F(ab') ₂
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.