

Product Datasheet

KAPPA LIGHT CHAIN antibody (orb436536)

Description	Goat polyclonal antibody to KAPPA LIGHT CHAIN
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human kappa light chains.
Target	KAPPA LIGHT CHAIN
Preservatives	None present. Borate buffered saline
Form/Appearance	F(ab)2 fragments of purified IgG - liquid
Concentration	0.5 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.