

Product Datasheet

KAPPA LIGHT CHAIN antibody (Biotin) (orb436537)

Description	Goat polyclonal antibody against KAPPA LIGHT CHAIN conjugated to Biotin
Species/Host	Goat
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Human kappa light chains.
Target	KAPPA LIGHT CHAIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	F(ab)2 fragments of purified IgG conjugated to Biotin - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.