

Product Datasheet

LAMBDA LIGHT CHAIN antibody (RPE) (orb436547)

Description	Goat polyclonal antibody against LAMBDA LIGHT CHAIN conjugated to RPE
Species/Host	Goat
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Human lambda light chains.
Target	LAMBDA LIGHT CHAIN
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	F(ab') ₂ fragments of purified IgG conjugated to R. Phycoerythrin - liquid
Concentration	0.25mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.