

Product Datasheet

IgG1 antibody (RPE) (orb436554)

Description	Goat polyclonal antibody against IgG1 conjugated to RPE
Species/Host	Goat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Mouse IgG1 paraproteins.
Target	IgG1
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.