

Product Datasheet

IgG2a antibody (AP) (orb436556)

Description	Goat polyclonal antibody against IgG2a conjugated to AP
Species/Host	Goat
Reactivity	Mouse
Conjugation	AP
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Immunogen	IgG2a paraproteins from BALB/c mice
Target	IgG2a
Preservatives	0.1% Sodium Azide (NaN ₃)50% Glycerol. TRIS buffered saline, 1mM MgCl ₂
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.