

Product Datasheet

IgG2a antibody (RPE) (orb436559)

Description	Goat polyclonal antibody against IgG2a conjugated to RPE
Species/Host	Goat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	IgG2a paraproteins from BALB/c mice
Target	IgG2a
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Concentration	0.5mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.