

Product Datasheet

IgG2b antibody (AP) (orb436561)

Description	Goat polyclonal antibody against IgG2b conjugated to AP
Species/Host	Goat
Reactivity	Mouse
Conjugation	AP
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Immunogen	Mouse IgG2b paraproteins.
Target	IgG2b
Preservatives	0.1% Sodium Azide (NaN ₃)50% Glycerol. 50mM TRIS1mM MgCl ₂
Form/Appearance	Purified IgG conjugated to alkaline phosphatase - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.