

Product Datasheet

IgG2c antibody (AP) (orb436565)

Description	Goat polyclonal antibody against IgG2c conjugated to AP
Species/Host	Goat
Reactivity	Mouse
Conjugation	AP
Tested Applications	ELISA
Immunogen	Native mouse IgG2c.
Target	IgG2c
Preservatives	0.09% Sodium Azide 0.2% Bovine Serum Albumin. 50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.