

Product Datasheet

IgG3 antibody (FITC) (orb436569)

Description	Goat polyclonal antibody against IgG3 conjugated to FITC
Species/Host	Goat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	Mouse IgG3 paraproteins.
Target	IgG3
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.