

Product Datasheet

IgM antibody (FITC) (orb436588)

Description	Goat polyclonal antibody against IgM conjugated to FITC
Species/Host	Goat
Reactivity	Human
Conjugation	FITC
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	Human IgM paraproteins.
Target	IgM
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.