

## Product Datasheet

### IgM antibody (RPE) (orb436592)

|                            |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody against IgM conjugated to RPE                                      |
| <b>Species/Host</b>        | Goat                                                                                        |
| <b>Reactivity</b>          | Human                                                                                       |
| <b>Conjugation</b>         | RPE                                                                                         |
| <b>Tested Applications</b> | FC                                                                                          |
| <b>Immunogen</b>           | Human IgM.                                                                                  |
| <b>Target</b>              | IgM                                                                                         |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> )Stabilizing agent (sucrose). Phosphate buffered saline |
| <b>Form/Appearance</b>     | F(ab') <sub>2</sub> fragment of IgG conjugated to R. Phycoerythrin (RPE) - liquid           |
| <b>Concentration</b>       | 0.25mg/ml                                                                                   |
| <b>Storage</b>             | +4°C DO NOT FREEZE                                                                          |
| <b>Note</b>                | For research use only                                                                       |
| <b>Isotype</b>             | Polyclonal IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                             |