

Product Datasheet

IgG antibody (Biotin) (orb436601)

Description	Goat polyclonal antibody against IgG conjugated to Biotin
Species/Host	Goat
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Normal human IgG.
Target	IgG
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(Fab')
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.