

Product Datasheet

KAPPA LIGHT CHAIN antibody (orb436605)

Description	Goat polyclonal antibody to KAPPA LIGHT CHAIN
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Human myeloma proteins with kappa light chains.
Target	KAPPA LIGHT CHAIN
Preservatives	None present. Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.