

Product Datasheet

IgG antibody (orb436611)

Description	Goat polyclonal antibody to IgG
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Normal rabbit IgG.
Target	IgG
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(Fab') (H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.