

Product Datasheet

IgG antibody (FITC) (orb436616)

Description	Rabbit polyclonal antibody against IgG conjugated to FITC
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Target	IgG
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	F(ab') ₂ fragment of IgG conjugated to Fluorescein Isothiocyanate Isomer I (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.