

Product Datasheet

IgG1 antibody (orb436621)

Description	Goat polyclonal antibody to IgG1
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Target	IgG1
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	F(ab') ₂ fragments of purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	HEAVY CHAIN
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.