

Product Datasheet

IgG antibody (HRP) (orb436630)

Description	Donkey polyclonal antibody against IgG conjugated to HRP
Species/Host	Donkey
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Purified goat IgG.
Target	IgG
Preservatives	0.05% Pro-clean 400.. Phosphate buffered saline.
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.