

Product Datasheet

IgG antibody (HRP) (orb436633)

Description	Rabbit polyclonal antibody against IgG conjugated to HRP
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, IH, IH-Fr, IHC-P
Immunogen	Purified rat immunoglobulin G.
Target	IgG
Preservatives	0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.