

Product Datasheet

IgG antibody (orb436634)

Description	Rabbit polyclonal antibody to IgG
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Mouse IgG.
Target	IgG
Preservatives	0.09% Sodium Azide. Borate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.