

Product Datasheet

IgG antibody (FITC) (orb436635)

Description	Sheep polyclonal antibody against IgG conjugated to FITC
Species/Host	Sheep
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC, IH-Fr
Immunogen	Rabbit IgG.
Target	IgG
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.