

Product Datasheet

APAAP COMPLEX antibody (orb436639)

Description	Mouse monoclonal antibody to APAAP COMPLEX
Species/Host	Mouse
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	IH, IH-Fr, IHC-P
Immunogen	Bovine alkaline phosphatase.
Target	APAAP COMPLEX
Preservatives	0.09% Sodium Azide. RPMI1640 medium
Form/Appearance	Complexes of monoclonal anti-alkaline phosphatase antibody with alkaline phosphatase
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Complex
Expiration Date	12 months from date of receipt.