

Product Datasheet

IgG antibody (orb436645)

Description	Goat polyclonal antibody to IgG
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Mouse IgG.
Target	IgG
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(RAT ADSORBED)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.