

Product Datasheet

IgG antibody (Biotin) (orb436648)

Description	Goat polyclonal antibody against IgG conjugated to Biotin
Species/Host	Goat
Reactivity	Hamster
Conjugation	Biotin
Tested Applications	ELISA, FC
Immunogen	Hamster IgG.
Target	IgG
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.8 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.