

Product Datasheet

IgG antibody (HRP) (orb436650)

Description	Goat polyclonal antibody against IgG conjugated to HRP
Species/Host	Goat
Reactivity	Hamster
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Hamster IgG.
Target	IgG
Preservatives	0.01%Thiomersal1%Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	0.8 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.