



Product Datasheet

Cytokeratin 19 antibody (orb43666)


Description	Mouse Monoclonal to Cytokeratin 19.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	Mammary organoids.
Target	Cytokeratin 19
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.Immunohistochemistry: Recommended dilution: 5-10 µg/ml.Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P08727
Entrez	3880
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.Immunohistochemistry: Recommended dilution: 5-10 µg/ml.Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of MCF-7 cells (red-filled) f...



Immunohistochemistry staining of human c...



Immunohistochemistry staining (paraffin ...