

Product Datasheet

IgG antibody (FITC) (orb436663)

Description	Goat polyclonal antibody against IgG conjugated to FITC
Species/Host	Goat
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, IH-Fr
Immunogen	Human IgG.
Target	IgG
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	F(ab') ₂ fragment of IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.7 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.