

Product Datasheet

IgG antibody (HRP) (orb436664)

Description	Goat polyclonal antibody against IgG conjugated to HRP
Species/Host	Goat
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, IH-Fr
Immunogen	Human IgG.
Target	IgG
Preservatives	0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	F(ab') ₂ fragments of purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	0.6 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.