

Product Datasheet

IgG antibody (RPE) (orb436665)

Description	Goat polyclonal antibody against IgG conjugated to RPE
Species/Host	Goat
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Human IgG.
Target	IgG
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	F(ab') ₂ fragment of IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.