

Product Datasheet

IgG2b antibody (orb436674)

Description	Rat monoclonal antibody to IgG2b
Reactivity	Rat
Tested Applications	ELISA
Target	IgG2b
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified myeloma protein - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	KAPPA
Isotype	IgG2b
Purity	Purified
Expiration Date	12 months from date of receipt.