

Product Datasheet

IgG1 antibody (orb436676)

Description	Mouse monoclonal antibody to IgG1
Reactivity	Mouse
Tested Applications	ELISA
Target	IgG1
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified monoclonal IgG1 - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Purity	Purified
Expiration Date	12 months from date of receipt.