

Product Datasheet

IgM antibody (orb436678)

Description	Pig monoclonal antibody to IgM
Reactivity	Porcine
Tested Applications	ELISA
Target	IgM
Preservatives	0.1% Sodium Azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified porcine immunoglobulin from serum - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgM
Purity	Purified
Expiration Date	12 months from date of receipt.