

Product Datasheet

IgG antibody (orb436685)

Description	Rabbit polyclonal antibody to IgG
Reactivity	Rabbit
Tested Applications	ELISA
Target	IgG
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified polyclonal rabbit IgG - liquid
Concentration	2 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Purity	Purified
Expiration Date	12 months from date of receipt.