


Product Datasheet

HLA-A2 antibody (orb43676)

Description	Mouse Monoclonal to HLA-A2.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	HLA-A2 solubilised by papain
Target	HLA-A2
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...