

## Product Datasheet

### Hrp-Stabilizing Diluent (orb436781)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Hrp-Stabilizing Diluent         |
| <b>Tested Applications</b> | ELISA, IH-Fr, IHC-P, WB         |
| <b>Target</b>              | HRP-STABILIZING DILUENT         |
| <b>Preservatives</b>       | 15ppm CMIT/MIT. Ready To Use    |
| <b>Form/Appearance</b>     | Ready to use - liquid           |
| <b>Storage</b>             | +4°C DO NOT FREEZE              |
| <b>Note</b>                | For research use only           |
| <b>Hazard Information</b>  | Non-Toxic                       |
| <b>Expiration Date</b>     | 12 months from date of receipt. |