

Product Datasheet

Antigen Retrieval Buffer, Universal (orb436807)

Description	Antigen Retrieval Buffer, Universal
Tested Applications	IHC-P
Target	ANTIGEN RETRIEVAL BUFFER, Universal
Preservatives	Bio-Rad's Antigen unmasking fluid is supplied as 3x concentrated fluid. Crystals in the fluid will dissolve during heating before use.
Storage	+25°C
Note	For research use only
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.