

Product Datasheet

ELISA PLASMA SAMPLE DILUENT (orb436810)

Description	ELISA PLASMA SAMPLE DILUENT
Tested Applications	ELISA
Target	ELISA PLASMA SAMPLE DILUENT
Preservatives	0.17% Proclin 30010% Bovine Serum Albumin5% EDTA. Phosphate buffered saline.
Form/Appearance	Ready to use - liquid.
Storage	+4°C DO NOT FREEZE
Note	For research use only
Expiration Date	12 months from date of receipt.