

Product Datasheet

Human ADIPOR2 ELISA Kit (orb436823)

Description	Human Adiponectin Receptor 2 (ADIPOR2) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q86V24
Sensitivity	0.24ng/mL
Entrez	79602
Expiration Date	Please enquire.