

Product Datasheet

Bovine ADP ELISA Kit (orb436825)

Description	Bovine Adiponectin (ADP) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q3Y5Z3
Sensitivity	32ng/mL
Entrez	282865
Expiration Date	Please enquire.