

Product Datasheet

Chicken ADP ELISA Kit (orb436827)

Description	Chicken Adiponectin (ADP) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6QWE7
Sensitivity	0.292ng/mL
Entrez	404536
Expiration Date	Please enquire.