

Product Datasheet

Rat ADP ELISA Kit (orb436831)

Description	Rat Adiponectin (ADP) Sandwich ELISA Kit
Reactivity	Rat
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8K3R4
Sensitivity	1.36ng/mL
Entrez	246253
Expiration Date	Please enquire.