

Product Datasheet

Rat AGTR2 ELISA Kit (orb436838)

Description	Rat Angiotensin II Receptor 2 (AGTR2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P35351
Sensitivity	0.051ng/mL
Entrez	24182
Expiration Date	Please enquire.