

Product Datasheet

Canine AngI ELISA Kit (orb436848)

Description	Canine Angiotensin I (AngI) competitive inhibition ELISA Kit
Reactivity	Canine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Sensitivity	9.28pg/mL
Expiration Date	Please enquire.