

Product Datasheet

Bovine AngII ELISA Kit (orb436850)

Description	Bovine Angiotensin II (AngII) competitive inhibition ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	8.46pg/mL
Expiration Date	Please enquire.