

Product Datasheet

Rat AngII ELISA Kit (orb436854)

Description	Rat Angiotensin II (AngII) competitive inhibition ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	24.31pg/mL
Expiration Date	Please enquire.