

Product Datasheet

Rat Bcl2L2 ELISA Kit (orb436864)

Description	Rat B-Cell CLL/Lymphoma 2 Like Protein 2 (Bcl2L2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O88996
Sensitivity	0.067ng/mL
Entrez	60434
Expiration Date	Please enquire.